

## List of consumables for chewing simulators

Picture	Part-No.	Description	Packing Unit / Qty
	0.02.011	Holder for Composite / Steatite sphere	50
	0.02.012	Socket for composite holder with specimen station connector	1
	0.01.011	Stainless steel antagonist tip radius 1 mm / 30° cone	8
	0.01.012	Stainless steel antagonist tip radius 2 mm / 30° cone	
	0.01.013	Stainless steel antagonist Ø8 mm sphere	8
	0.01.014	Stainless steel antagonist Ø6 mm sphere	
	0.01.015	Stainless steel antagonist Ø4 mm sphere	
	0.01.008	Steatite sphere Ø6 mm	50
	0.02.013	Multipurpose antagonist (e.g. for natural teeth) inner Ø15 mm, height 11 mm	8
	0.01.001	Antagonist adapter Ø16 – Ø5 mm	1

Picture	Part-No.	Description	Packing Unit / Qty
	0.01.010	Antagonist adapter Ø16 – M10 thread To use specimen container as antagonist	1
	0.02.006-K	Specimen container, inner Ø30 mm with metal pin and specimen station connector	1
	0.02.015-K	Specimen container, inner Ø42 mm with metal pin and specimen station connector	1
	0.02.015-K	Specimen container, inner Ø76 mm with specimen station connector and set screws	1
	0.02.005	Specimen station connector	1
	0.02.016-K	Specimen holder for DIN 14801 testing, 30° load angle with fixation screw and specimen station connector	1
<b>Single-use consumables</b>			
	0.01.016	ABS-holder with steatite sphere Ø6 mm	1
	0.01.017	ABS-holder with carbide sphere Ø4 mm	8
	0.01.018	ABS-holder with carbide sphere Ø6 mm	
	0.01.019	ABS-holder with carbide sphere Ø8 mm	
	0.02.019	ABS-specimen container, inner Ø25 mm	1
	0.02.020	ABS-specimen container, inner Ø30 mm	
	0.02.021	ABS-specimen container, inner Ø42 mm	
	0.02.022	ABS-specimen container, inner Ø76 mm	
	0.02.023	Wing screw for single-use articles	1

**Remark:** Carbide means tungsten carbide